

Features

- Modular configuration
- Wall mount: One to six (1-6) 1 liter bag capacity
- Countertop: One to three (1-3) 1 liter bag capacity
- Pole mount: One (1) 1 liter bag capacity
- Anodized aluminum heat plate
- Lightweight, medical grade plastic housing
- Intelligent bag sensors
- Dual temperature sensors
- Heats to 40°C (104°F) in 30 minutes or less
- Accuracy of +0°/-2°C (+0°/-3°F)



ENTHERMICS Medical Systems

W164 N9221 Water Street
P.O. Box 443
Menomonee Falls, Wisconsin
53052-0443 U.S.A.

Phone: (262) 251-8356
(800) **TO-B-WARM** U.S.A./CANADA
(800) 862-9276 U.S.A./CANADA

Fax: (262) 251-7067
(800) 329-8744 U.S.A./CANADA

Email: generalinfo@enthermics.com

Web site: www.enthermics.com

ISO 13485:2003 Certified

#E123-06/10

Experience the
Entermics difference:
quality beyond comfort.



ivNow™

Injection/Irrigation Fluid Warmers
Patent WO 2010/030746



Modular Thermal Fluid Warming

ivNow modular fluid warmers put warmed fluid at your fingertips where you need it, when you need it. ivNow fluid warmers are small, convenient units that can be placed in a number of different combinations and environments. Countertop, wall mounted and mobile units are all available depending on your particular requirements.

ivNow quickly warms and maintains safe temperatures of intravenous and irrigation fluids. The specially contoured warming surface cradles one (1) liter solution bags. Each cavity is controlled by an individual control and L.E.D. display. A sensor in the heating plate detects the presence of a bag and engages the heating mechanism to quickly warm the fluid. Two temperature sensors continuously read the temperature of the fluid bag to engage the heater as necessary to maintain the fluid temperature within +0/-2°C (+0/-3°F) of the set point temperature.

ivNow warmers ensure your fluids are kept at the proper temperature and will always be immediately available at or close to body temperature rather than room temperature without the added expense and set-up time required with in-line warming systems.

See ivNow in action at www.youtube.com/user/Enthermics and www.facebook.com/pages/Enthermics

Available Configurations

ivNow-1 (one (1) 1 liter bag capacity) can be placed directly on a countertop, mounted on an optional countertop pedestal stand, mounted on a mobile equipment pole stand, or mounted on a wall using brackets.

ivNow-2 (two (2) 1 liter bag capacity) can be placed directly on a countertop, mounted on a countertop pedestal stand, or mounted on a wall using brackets.

ivNow-3 (three (3) 1 liter bag capacity) must be mounted on a countertop pedestal stand or mounted on a wall using mounting brackets.

ivNow-4 (four (4) 1 liter bag capacity) must be mounted on a wall using mounting brackets.

ivNow-5 (five (5) 1 liter bag capacity) must be mounted on a wall using mounting brackets.

ivNow-6 (six (6) 1 liter bag capacity) must be mounted on a wall using mounting brackets.

Note: Individual ivNow units cannot be stacked in the field. Multiple cavity units are available from the factory in the configurations above.

To learn more, contact your Entermics representative at (800) TO-B-WARM or visit us at www.enthermics.com.

