

# Level 356 Processor



MD Version

**No Plumbing Required!**  
**No Darkroom Required with Daylight Loader Option!**  
**Ultra-Compact Design**  
**Takes up to 35cm Wide Films**  
**Visible Chemical Level Tubes**  
**Front Fill & Front Drain**  
**Front Temperature Control**  
**Energy-Saving Standby Switch**  
**Uni-Rack Simplifies Cleaning**  
**Convenient Front Panel Fuse**  
**1 Year Warranty**



MD-D Version

## SIMPLE INSTALLATION, OPERATION AND MAINTENANCE

### EASY OPERATION



Fill and drain the processor using the level tubes in the front of the processor.

The Temperature dial and Run/Standby dial are easily accessible from the front of the processor.

### EASY DAILY MAINTENANCE



Uni-Rack makes for easy daily roller maintenance.

## SPECIFICATIONS

Film Size	30mm X 30mm Minimum, Up to 356mm Wide X Any Length
Process Time	Fixed 4.5 Minutes
Process Temperature	Adjustable from 20°C to 40°C
Tank Capacity	1.5 liters Each of Developer, Fixer, and Wash
Circulation	Independent Pumps for Developer, Fixer, and Wash
Line Voltage	115VAC or 230VAC, 8 Amps
Dimensions	550mm Wide X 400mm Deep X 540mm Tall with Daylight Loader 550mm Wide X 400mm Deep X 350mm Tall without Daylight Loader
Net Weight (no fluids)	Approximately 20 kilograms.
Shipping Size & Weight	610mm X 420mm X 620mm, 23 kilograms

**J.D. Honigberg International, Inc.**

85 Revere Drive, Suite D  
 Northbrook, IL 60062 USA  
 (847) 412-0200 Fax: (847) 412-0034

**Visit our website at**  
**[www.jdhmedical.com](http://www.jdhmedical.com)**  
**or email us at**  
**[medical@jdhintl.com](mailto:medical@jdhintl.com)**