

# CASTING & SPLINTING ALUMINUM / STEEL CARTS

Stainless Steel Casting and Splinting System with pump pole bracket, bucket holder with vinyl mat, aluminum towel bar, folding shelf, swivel tourniquet/pump pole, top shelving system, side storage cabinet, scissors cable and adjustable drawer dividers



**6025-TC**



**6035TC**

Painted Casting and Splinting System with steel cabinet and aluminum drawers. Includes a pump pole bracket, bucket holder with vinyl mat, aluminum towel bar, folding shelf, swivel tourniquet/pump pole, top shelving system, side storage cabinet, scissors cable and adjustable drawer dividers



Stainless Steel Casting and Splinting System with pump pole bracket, bucket holder with vinyl mat, aluminum towel bar, folding shelf, swivel tourniquet/pump pole, side storage cabinet, scissors cable and adjustable drawer dividers



**6025**



**6035**

Painted Casting and Splinting System with steel cabinet and aluminum drawers. Includes a pump pole bracket, bucket holder with vinyl mat, aluminum towel bar, folding shelf, swivel tourniquet/pump pole, side storage cabinet, scissors cable and adjustable drawer dividers

Stainless Steel Casting Cart with pump pole bracket, bucket holder with vinyl mat, aluminum towel bar and adjustable drawer dividers



**6020**



**6030**

Painted Casting Cart with steel cabinet and aluminum drawers. Includes a pump pole bracket, bucket holder with vinyl mat, aluminum towel bar and adjustable drawer dividers

## Overall Dimensions

43" H x 44" W x 21"D  
(109 cm x 112 cm x 53 cm)

## Shipping Dimensions

45" H x 55" W x 29" D  
(114 cm x 140 cm x 74 cm)

## Internal Drawer Dimensions

16.5" D x 23.5" W  
(42 cm x 59.7 cm)

## Working Surface on top of cart

41.5" H (105.5 cm) above ground

## Net Weight

**6025-TC** 219 lbs (100 kg)  
**6035-TC** 219 lbs (100 kg)  
**6025** 209 lbs (95 kg)  
**6035** 209 lbs (95 kg)  
**6020** 182 lbs (83 kg)  
**6030** 182 lbs (83 kg)

## Shipping Weight

28 lbs (127 kg)  
 28 lbs (127 kg)  
 270 lbs (122 kg)  
 270 lbs (122 kg)  
 240 lbs (109 kg)  
 240 lbs (109 kg)

## Drawers:

- 1 ea 3" top drawer with black vinyl mat
- 5 ea 3" drawer with 4 adjustable dividers per drawer
- 2 ea 6" drawers with 4 adjustable dividers per drawer

**Standard Features include:**

- All Stainless steel construction or Painted Steel Cabinet and Base
- Stainless steel bucket holder with black vinyl mat (holds two casting buckets)
- Aluminum towel bar on the back of the cart
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) premium rubber single-wheel caster set, full swivel – one with brake, one directional
- 24-inch (61 cm) stainless steel ruler, mounted on the top of the cart
- Heavy-duty ball-bearing drawer slides
- Cart is delivered fully assembled with 5 Year Warranty
- Aluminum Drawers (6035 / 6030)
- Certified “Latex Free” (6035 / 6030)
- Hard-baked powder coat finish in 12 colors (6035-TC / 6035 / 6030)

**Standard Features ONLY for carts model 6025-TC, 6035-TC, 6035, 6030**

- Folding shelf (6035-TC / 6035 / 6030)
- Side storage cabinet (6035-TC / 6035 / 6030)
- Swivel tourniquet / pump pole (6035-TC / 6035 / 6030)
- Scissors cable (6035-TC / 6035 / 6030)

**Standard Features ONLY for carts model 6025-TC, 6035-TC**

- Top compartment for splinting supplies (602106)

**Thirteen available color choices:**

- Hammer Beige (HB), Gloss Beige (BG), Burgundy (BU), Light Blue (LB), Navy (NY), Light Gray (GY), Red (RD), Teal (TL), Mauve (MV), Yellow (YL), Hammer Tone Navy (HTDB), Hammer Tone Gray (HTGY), Hammer Tone Green (HTGN)

**Optional Accessories for AL16CS, AL20CS, AL24CS, AL28CS:**

- 602001 Stainless Steel Drop Shelf, Cast Cart
- 602002 Stainless Steel Side Storage Cabinet, Cast Cart
- 602003 Stainless Steel Bucket Bracket With Swivel Tourniquet Pump Pole
- 602003R Retrofit Bucket Bracket With Swivel Tourniquet Pump Pole, For Field Installation
- 602004 Scissors Security Cable, Cast Cart
- 602005 Clear Label Holders, Cast Cart, Pk-10
- 602006 Top Compartment/Organizer, Cast Cart - Holds 3M Supplies
- 602106 Top Compartment/Organizer, Cast Cart - NEW - Holds BSN and 3M Supplies
- 602007 Swivel Tourniquet Pump Pole
- 602012 Black Vinyl Drawer Mat, Cast Cart
- 602013 3" (7.6 cm) Drawer, Stainless Steel, Cast Cart
- 602014 3" (7.6 cm) Drawer Divider Kit, Cast Cart
- 602016 6" (15.2 cm) Drawer, Stainless Steel, Cast Cart
- 602017 6" (15.2 cm) Drawer Divider Kit, Cast Cart