

KEY LOCKING STEEL PUNCH CARD MEDICATION CARTS



VL280PC
280 Card capacity



SL480PC
480 Card capacity



SL600PC & SL600PCSP (Special Package)
600 Card capacity

| | Overall Dimensions | Internal Drawer Dimensions | Net Weight | Shipping Weight |
|-----------------|--|---|---|------------------------|
| VL280PC: | 44.75"H x 29.5" W x 22.5"D (114 cm x 75 cm x 57 cm) | 17" D x 23" W (43.2 cm x 58.4 cm) | 140 lbs (64 kg) | 175 lbs (80 kg) |
| SL480PC | 45.75"H x 36" W x 23.5"D (116 cm x 91 cm x 60 cm) | <u>Left Side:</u> 17"D x 20.5"W (43.2 x 52 cm) <u>Right Side:</u> 17"D x 6.62"W (43.2 x 16.8 cm) | 170 lbs (77 kg) | 215 lbs (98 kg) |
| SL600PC | 45.75"H x 43" W x 23.5"D (116 cm x 109 cm x 60 cm) | <u>Left Side:</u> 17"D x 20.5"W (43.2 x 52 cm) <u>Right Side:</u> 17"D x 6.62"W (43.2 x 16.8 cm) | 184 lbs (84 kg) | 229 lbs (104 kg) |
| | Shipping Dimensions | | Working Surface on top of cart | |
| VL280PC | 48" H x 25" W x 33" D (122 cm x 63.5 cm x 84 cm) | | 43" H (109 cm) above ground | |
| SL480PC | 52" H x 38" W x 26" D / (132 cm x 97 cm x 66 cm) | SL600PC 52" H x 46" W x 26" D (132 cm x 117 cm x 66 cm) | 44" H (112 cm) above ground (SL480/600PC) | |

Standard Features include:

- Pull out Side Shelf (VL280PC)
- Bumper Frame with extended stabilizing front base, 4 Corner Bumpers (VL280PC) or wrap around bumper (SL480/600PC)
- **VL280PC:** Five-inch (12.7 cm) premium rubber single-wheel casters, full swivel – one with brake, one directional
- **SL480PC & SL600PC:** 5" (12.7 cm) sealed ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional, easy to maneuver over all floor types; tile, carpet, linoleum, wood and concrete
- Painted steel top rail/push handle on VL280PC
- Removable plastic top for VL280PC or plastic replaceable top with integrated push handles for SL480PC & SL600PC)
- Full extension Ball-Bearing Drawer Slides (self closing for SL480PC & SL600PC)
- Cart is delivered fully assembled with 5 years Warranty
- Certified "Latex Free"
- Hard-baked powder coat finish

Features only on SL480PC & SL600PC models

- Dual Column Medication Delivery Cart with 480 / 600 Punch Card capacity
- Dual Locking System complies with regulatory security standard
- The "slam action" lock on the Narcotics drawer allows the drawer to be closed without using the key.
- An internal drawer cover keeps the narcotics drawer contents secure, even if all other drawers are removed

Drawers:

- VL280PC:**
- 1 x 3-inch (7.6 cm) Multipurpose storage drawer
 - 2 x 9-inch (22.9 cm) Punch card drawers with two adjustable dividers
 - 1 x 9-inch (22.9 cm) drawer with locking narcotics compartment and one adjustable divider

SL480/600PC: Right Column Drawers:

- 1 each – 4 inch (10.1 cm) Multipurpose storage drawer
- 1 each – 9.75 inch (24.8 cm) Locking Narcotics Drawer with 40 cards capacity
- 2 each – 9.75 inch (24.8 cm) Multipurpose Storage Drawer with 40 cards capacity per drawer

Left Column Drawers:

- 1 each – 4 inch (10.1 cm) Multipurpose storage drawer
- 3 each – 9.75 inch (24.8 cm) Punch card drawers with 160 cards capacity per drawer

Thirteen available color choices:

- Hammer Beige (HB), Gloss Beige (BG), Burgundy (BU), Light Blue (LB), Navy (NY), Light Gray (GY), Red (RD), Teal (TL), Mauve (MV), Yellow (YL), Hammer Tone Navy (HTDB), Hammer Tone Gray (HTGY), Hammer Tone Green (HTGN)

Optional Accessories for VL280PC, SL480PC, SL600PC:

- 683802 Sharps Container And Mounting Bracket, 3 gallon
- 40636 Sharps Container Only, 3 gallon
- 40638 Sharps Bracket for 3 gal Sharps Container
- 40639 Locking Sharps Bracket for 3 gal Sharps Container
- 683812 Sharps Container With Locking Bracket, 3 gallon
- 680440 Locking Sharps Holder, Includes 5 Qt. Sharps Container (40637)
- 680445 Locking Sharps Holder with Glove Box Holder (40633), Includes 5 Qt. Sharps Container (40637)
- 40633 NEW Glove Box holder only for 680440/40637 Sharps Holder
- 40635 Locking Sharps Box only, 5 Qt capacity (part of 680440/680445 accessory)
- 40637 Sharps Container only, 5 Qt, For 680440 Holder
- 681702 Sharps Bracket, Universal
- 68161 Aluminum Push Handle, 16"
- 680416 Chart Holder, Steel
- 680418 Cup Dispenser - 2 To 5 Oz Cups
- 680401 IV Pole - 2 Prong
- 40106 Top, Plastic Removable, Classic Series
- 680415 Tape/Label Dispenser
- 68160 Top Rail/Push Handle, Stainless Steel
- 68160SPL Top Rail/Push Handle, Split, Stainless Steel
- 68161EX Top Rail Subst, Replace 68161 Push Handle W/68160
- 680407 Utility Basket, Large, With Rail Clips
- 680447 Utility Basket, Large, With Mounting Brackets
- 40257 Utility Basket, Large
- 40256 Utility Basket, Small
- 680431 Utility Basket, Small, With Nylon Mounting Clamps
- 680402 Utility Hooks (Set Of 4)
- 680428 Waste Container Without Cover and mounting bracket, Plastic, 3 Gal.
- 40717 Waste Container only, Plastic
- 680408 Waste Container Without Cover, Plastic, 2 gallon
- 680409 Waste Container With Cover, Steel
- 680410 Waste Container Without Cover, Steel
- 640360 Cup and Straw Organizer, White Styrene Plastic
- 681501 3" Economy Exchange Tray, Each
- 681502 6" Economy Exchange Tray, Each
- 680506 6/9/12" Drawer Divider Set with Three Adjustable dividers, factory installed only
- 40582 M/S Cart Drawer Divider Tray Only
- 68530-P1 M/S Cart Drawer Divider Tray Pack 1 with 3 Long, 5 Short. 2 (1/2)L, 2 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- 68530-P2 M/S Cart Drawer Divider Tray Pack 2 with 5 Long, 7 Short. 4 (1/2)L, 4 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- 680301 Side Mounted Drop Shelf With Lip
- 680302 Side Mounted Drop Shelf Without Lip
- 68175 UPGRADE 9" Punch Card Drawer to Add Locking Narcotics Box
- 40616 Patient Card Divider, Each
- 640616 Patient Card Divider, Pack Of 25
- 68640-P1 4w Drawer Divider Tray Pack 1 with 3 Long, 5 Short. 2 (1/2)L, 2 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- 68640-P2 4w Drawer Divider Tray Pack 2 with 5 Long, 7 Short. 4 (1/2)L, 4 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- 68630-P1 3w Drawer Divider Tray Pack 1 with 3 Long, 5 Short. 2 (1/2)L, 2 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- 68630-P2 3w Drawer Divider Tray Pack 2 with 5 Long, 7 Short. 4 (1/2)L, 4 (1/2)S, 1-2", 1-3", 1-4", 1-5"
- 680401MC IV Pole For Medication Cart
- **VL280PCSP** Specialty Package with Waste Container, Locking Sharps, Adjustable Dividers in top drawer Tray
- **SL480PCSP** Special Pack with Waste Container, Locking Sharps, Pullout shelf, Adjustable Dividers in top drawer Tray.
- **SL600PCSP** Specialty Package with Waste Container, Locking Sharps, Adjustable Dividers in top drawer Tray
- **VL280MPC** Same as model VL280PC but with keyless push button entry system instead of key lock entry system
- **VL280MPCSP** Same as model VL280PCSP but with keyless push button entry system instead of key lock entry system.

