



MEDICATION CASSETTE CARTS OR MEDICATION BIN CARTS?

PATIENT BIN CARTS

Harloff's Bin carts have the patient bin drawers in the cart, resting directly inside the drawers (see picture below) and each bin is pulled out individually. Locking system is via key lock but some models can have an optional push button lock system instead.



The Bin carts are available in various sized cart frames and come in 4 versions with different bin sizes and quantity.

- 1) Model VL12BIN5: 12 individual bins cart (each bin is 5 inches / 12.7 cm wide)
- 2) Model VL18BIN: 18 individual bins car (each bin is 3.5 inches / 8.9 cm wide)
- 3) Model VL20BIN5: 20 individual bins cart (each bin is 5 inches / 12.7cm wide)
- 4) Model VL24BIN5: 24 individual bins cart (each bin is 5 inches / 12.7cm wide); optional push button lock.
- 5) Model VL30BIN3: 30 individual bins cart (each bin is 3.5 inches / 8.9 cm wide)
- 6) Model VL36BIN3: 36 individual bins cart (each bin is 3.5 inches / 8.9 cm wide); optional push button lock.
- 7) Model SL28BIN5: 28 individual bins cart (each bin is 5 inches / 12.7 cm wide)
- 8) Model SL36BIN5: 36 individual bins cart (each bin is 5 inches / 12.7 cm wide)
- 9) Model SL42BIN3: 42 individual bins cart (each bin is 3.5 inches / 8.9 cm wide)
- 10) Model SL54BIN3: 54 individual bins cart (each bin is 3.5 inches / 8.9 cm wide)



Bin drawer with the patient bins



Bin dividers inside patient bins



VL Bin cart series



SL Bin cart series

CASSETTE CART SYSTEM

Harloff's Cassette cart system consists in a cassette holding inside it one or two layers of Bins. This system allows for easy exchange, where multiple bin drawers can be removed from the cart at once with one hand (see picture below).

The cassette with 1 layer of bins holds a total of 4 bins. The cassette with 2 layers of bins holds a total of 8 bins (see picture below)

The cassette carts are only available in one sized cart frame and can come in 4 versions for the total number of bins. There is only 1 size for the bins (4.2 inches / 10.7 cm wide).

- 1) Model AL16CS: 16 bins cart (2 cassettes of 8 bins each)
- 2) Model AL20CS: 20 bins cart (2 cassettes of 8 bins each & 1 cassette of 4 bins)
- 3) Model AL24CS: 24 bins cart (3 cassettes of 8 bins each)
- 4) Model AL28CS: 28 bins cart (3 cassettes of 8 bins each & 1 cassette of 4 bins)



1 & 2 bin layer cassettes



Cassette Cart series



Cassette being taken out of cart drawer