

# MRI SAFE EMERGENCY & ANESTHESIA ALUMINUM CARTS



**MR6B (Emergency)**



**MR7B (Emergency)**



**MR6K (Anesthesia)**



**MR7K (Anesthesia)**

Overall Dimensions	Internal Drawer Dimensions	Net Weight		Shipping Weight
45.75" H x 30.5" W x 23.75"D (116 cm x 77.5 cm x 60 cm)	17" D x 23" W (43.2 cm x 58.4 cm)	<b>MR6B</b>	126 lbs (57 kg)	156 lbs (71 kg)
		<b>MR7B</b>	126 lbs (57 kg)	156 lbs (71 kg)
<b>Shipping Dimensions</b> 52" H x 26" W x 33" D (132 cm x 66 cm x 84 cm)	<b>Working Surface on top of cart</b> 44" H (112 cm) above ground	<b>MR6K</b>	126 lbs (57 kg)	156 lbs (71 kg)
		<b>MR7K</b>	126 lbs (57 kg)	156 lbs (71 kg)

### Standard Features include:

- Tall Cabinet with Breakaway locking system for emergency carts MR6B & MR7B
- Tall Cabinet with Key locking system for anesthesia carts MR6K & MR7K
- Vinyl padded full wrap around steel Bumper
- Five-inch (12.7 cm) premium ball bearing rubber single-wheel casters, full swivel – one with brake, one directional
- Stainless Steel Pull-Out Shelf
- Pre-threaded accessory mounting holes
- Replaceable Plastic Top with molded-in Push Handles
- Low Ferrous Ball-Bearing Drawer Slides
- Cart is delivered fully assembled with a 5 years Warranty
- Certified "Latex Free"
- Hard-baked powder coat finish

### Drawers:

#### Carts MR6B & MR6K

- 1 x 4 inch (10.2 cm) Multipurpose storage drawer
- 2 x 3.25 inch (8.3 cm) Multipurpose storage drawer
- 2 x 6.5 inch (16.5 cm) Multipurpose storage drawer
- 1 x 9.75 inch (24.8 cm) Multipurpose storage drawer

#### Carts MR7B & MR7K

- 1 x 4 inch (10.2 cm) Multipurpose storage drawer
- 4 x 3.25 inch (8.3 cm) Multipurpose storage drawer
- 1 x 6.5 inch (16.5 cm) Multipurpose storage drawer
- 1 x 9.75 inch (24.8 cm) Multipurpose storage drawer

### Thirteen available color choices:

- Hammer Beige (HB), Gloss Beige (BG), Burgundy (BU), Light Blue (LB), Navy (NY), Light Gray (GY), Red (RD), Teal (TL), Mauve (MV), Yellow (YL), Hammer Tone Navy (HTDB), Hammer Tone Gray (HTGY), Hammer Tone Green (HTGN)
- Shelf painted Gloss White with MRI Sticker

### Optional MR Safe Accessories for MR6B, MR7B, MR6K, MR7K:

- **MR6B-EMG** Cart MR6B with IV Pole, Drop Shelf, 4 Utility Hooks, Cardiac Board, O2 Tank holder, 2 drawer Trays with Dividers
- **MR7B-EMG** Cart MR7B with IV Pole, Drop Shelf, 4 Utility Hooks, Cardiac Board, O2 Tank holder, 2 drawer Trays with Dividers
- **MR6K-MAN** Cart MR6K with 3 Tier raised Back Rail system, 4 Side Rails, 2 Full Size drawer Trays with Dividers
- **MR7K-MAN** Cart MR7K with 3 Tier raised Back Rail system, 4 Side Rails, 2 Full Size drawer Trays with Dividers

**Optional MR Safe Accessories for MR6B, MR7B, MR6K, MR7K:**

- **MR-RBR3** 3 Tier Raised Back Rail for MRI compatible carts
- **MR-RBR2** 2 Tier Raised Back Rail for MRI compatible carts
- **MR-SR** Side Rail for MRI compatible carts
- **MR-CARDBRD** Cardiac Board & Mounting brackets for MRI compatible carts
- **MR-RLHK** Rail Hook for MRI compatible carts
- **MRIV2** IV Pole for MRI compatible carts
- **MR-O2** Brackets for O2 Tank for MRI compatible carts
- **MR-POSHELF** Replacement Pullout Shelf for MRI compatible carts
- **MR-UHOOKS4** Utility Hooks (set of 4) for MRI compatible carts
- **MR4TILT** 4 Compartment Tilt Bin with mounting clip for MRI compatible carts
- **MR-5TILT** 5 Compartment Tilt Bin with mounting clip for MRI compatible carts
- **MR-6TILT** 6 Compartment Tilt Bin with mounting clip for MRI compatible carts
- **MR-3GWASTE** Waste Container, 3 Gallon, with mounting bracket for MRI compatible carts
- **MR-8QWASTE** Waste Container, 8 Quart, with mounting bracket for MRI compatible carts