

# STEEL ANESTHESIA MALIGNANT HYPERTHERMIA CARTS



**6640**

**6641** (Specialty Package)

Overall Dimensions	Internal Drawer Dimensions	Net Weight	Shipping Weight
41.5" H x 32" W x 22"D (105 cm x 81 cm x 56 cm)	17" D x 23" W	<b>6640</b> 146 lbs (66.4 kg)	181 lbs (82.3 kg)
44.5" H x 32" W x 22"D (113 cm x 81 cm x 56 cm)	(43.2 cm x 58.4 cm)	<b>6641</b> 154 lbs (70 kg)	199 lbs (90.5 kg)

**Shipping Dimensions**  
49" H x 39" W x 28" D  
(124 cm x 99 cm x 71 cm)

**Working Surface on top of cart**  
40.5" H (103 cm) above ground

## Standard Features include:

- Tall Cabinet with breakaway lock seals (pack of 50)
- Separate breakaway locks on both the refrigerator and the two storage drawers
- 0.045 cubic meter refrigerator / freezer with hospital grade plug (120 Volts only for refrigerator and plug)
- 10 Glycol packs
- LED-illuminated digital thermometer with cart power safety feature
- Vinyl padded full wrap around steel Bumper
- Five-inch (12.7 cm) premium ball bearing rubber single-wheel casters, full swivel – one with brake, one directional
- Aluminum Push handle / side rail
- Removable plastic top
- Full extension Ball-Bearing Drawer Slides
- Cart is delivered fully assembled with a 5 years Warranty
- Certified "Latex Free"
- Hard-baked powder coat finish

## Specialty Package for 6641 include:

- IV pole
- Stainless steel Rail / Push handle
- Side mounted drop shelf
- Full drawer Tray organizer with adjustable dividers for top drawer

## Drawers:

- 1 x 3-inch (7.6 cm) Multipurpose storage drawer
- 1 x 6-inch (15.2 cm) Multipurpose storage drawer

## Thirteen available color choices:

- Hammer Beige (HB), Gloss Beige (BG), Burgundy (BU), Light Blue (LB), Navy (NY), Light Gray (GY), Red (RD), Teal (TL), Mauve (MV), Yellow (YL), Hammer Tone Navy (HTDB), Hammer Tone Gray (HTGY), Hammer Tone Green (HTGN)

## **Features and Benefits:**

- Some patients under anesthesia experience a sudden and dangerous rise in blood temperature, a condition called Malignant Hyperthermia (MH). Healthcare facilities that administer volatile inhalational agents known to trigger MH must have emergency equipment immediately available for the treatment of MH (Accreditation Association for Ambulatory Care). Surgery centers must have immediately available, in an organized cart, all items required for treatment of MH (American Society of Anesthesiologists and the American Association of Nurse Anesthetists).
- The Harloff Malignant Hyperthermia Cart will demonstrate JCAHO compliance and/or best practice for all OR suites, Level 1 Trauma Units, Surgery Centers, and Ambulatory Facilities.
- The Cart increases patient safety and departmental efficiency by providing a convenient, portable and economic solution for storage and transport of cooled glycol packs, cold saline IV solutions, and dantrolene (Dantrium) used to stabilize MH. Cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most likely needed.
- External digital thermometer is visible at a distance in all lighting condition, and allows continual monitoring of the supplies' temperature.
- Cart ships with wheels assembled, so the cart does not need to be inverted to install the wheels, thereby protecting the refrigeration components.