



750 STRETCHER



Pneumatic Assist Fowler

Pneumatic cylinder provides additional assist for the safe and easy raising and lowering of the patient.

Quick-Release Oxygen Tank Holder

Pedigo's Quick-Release O₂ Tank Holder is conveniently located under the litter reducing caregiver bending. (Patent Pending)

Dual Side Mounted Hydraulic Foot Pedals

Hydraulic Foot Pedals are located on each side of the stretcher within the patient care zone.



5th Wheel Steering

Optional 5th Wheel Steering increases maneuverability and control.

Quick-Release O₂ Holder

Ergonomically located for easy access in the replacement or loading of O₂ tanks. (Patent Pending)



750 STRETCHER



Pedigo's NEW 750 Series Stretcher allows for the transfer of patients easily and comfortably with a load capacity up to 750 pounds. Superior stability, mobility, and braking capabilities make the 750 ideal for ED, OR, PACU or any other department in your facility.



Tuck-Away Side Rails

Full length stainless steel side rails feature integrated push-pull handles with a latch release mounted at all four corners of the stretcher.

Hand Control Trendelenburg

Mounted at the foot end of the stretcher, the unique Hand Control Trendelenburg with pneumatic assist easily moves from horizontal to positive or reverse positions. Pedigo's unique mechanical memory lock perfectly returns the angle of the litter to the horizontal position.

Side Rail Latch Release

Located at all four corners of the stretcher, the convenient hand latch provides easy one-handed access to the release mechanism.



Brake & Steer Foot Pedals

Ergonomically designed, the Brake and Steer mechanism provides effortless heel and toe activation at all four corners.



STANDARD FEATURES

- 750 lbs. weight capacity
- 29.5" wide (750-N) or 33.5" wide (750-W)
- 23" low height
- Tuck-away side rails with integrated push/pull handles and latch release
- Quick-Release O₂ holder included (Patent Pending)
- Dual side-mounted hydraulic raise/lower foot pedals
- IV pole receptacles at every corner
- Four restraint strap positions
- Single hydraulic cylinder design
- Pneumatic assist positive and negative Trendelenburg
- Pneumatic assist Fowler backrest
- Brake and steer foot pedals, with heel and toe activation
- True Directional Steering or optional 5th Wheel Steering
- 8-inch easy mobility casters
- Mattress sold separately. See Options for sizes.

OPTIONS

- 5th Wheel Steering
- Knee Flex/Leg Lift
- Mattress Options:
 - 3" Standard Mattress
 - 3" Premium Mattress
 - 4" Standard Mattress
 - 4" Premium Mattress
 - 4" Pressure Care Mattress
- IV Pole Options:
 - Fixed
 - Adjustable / Removable
 - Adjustable / Attached
 - Fold Down
- Side Rail Pad
- Stirrups
- Stirrups Adapter
- Vertical Oxygen Tank Holder
- Monitor Tray Board
- Foot Board/Chart Holder
- Restraint Straps (Velcro)
- Restraint Straps (Buckle)

SPECIFICATIONS

Model Number	750-N	750-W
Overall Length	83" (210.8 cm)	83" (210.8 cm)
Overall Width	29.5" (74.9 cm)	33.5" (85 cm)
Weight Capacity	750 lbs. (340.2 kg)	750 lbs. (340.2 kg)
Height Range		
Low	23" (58.4 cm)	23" (58.4 cm)
High	32.75" (83.2 cm)	32.75" (83.2 cm)
Litter Positioning		
Fowler Backrest	0°-90°	0°-90°
Trendelenburg	+12°	+12°
Reverse Trendelenburg	-7°	-7°
Patient Surface	77" x 27" (195.6 x 68.6cm)	77" x 31" (195.6 x 78.7cm)
Side Rails	14" x 76" (35.6 x 193 cm)	14" x 76" (35.6 x 193 cm)
Casters	8" (20.3 cm)	8" (20.3 cm)
Warranty	1 year parts and labor	1 year parts and labor



ALL PEDIGO PRODUCTS ARE PROUDLY MADE IN THE USA