

The NEW 5240 Exam Table

same great quality; new great look



The new 5240 examination table maintains the tradition of strength and quality of the original UMF Signature Series, but with an updated look. • Reversible drawers eliminate the need of choosing a right or left handed table and provide flexibility of room layout. • The recessed paper storage area provides for extra storage and convenience. • Seamless molded drawers provide protection from messy spills. • Radiused corners give a new look to a traditional design. • The Ultra Comfort Seamless top is easy to clean and provides security and comfort for patients up to 500 pounds. • All other features of the original Signature Series table including fully-articulating back section, all steel footrest, easy glide comfort stirrups and quiet roller bearing drawer glides are incorporated in this new design. Of course, the new table comes with the same 5-year warranty as the original 5100 series.

*UMF Quality, Design
and Innovation*


EXPECT MORE.

See other side for standard features, options and dimensions.



The NEW 5240 Exam Table

same great quality; new great look



Reversible drawers eliminate the need of choosing a right or left handed table and provide flexibility of room layout.



Easy-glide comfort stirrups

Standard Features

- UMF Ultra Comfort Top (vacuum formed and seamless)
- Easy-glide comfort stirrups (3 positions)
- Upholstered removable leg rest
- Pull-out all steel non-slip step (7" high)
- Paper roll holder
- Plastic 9"x9" (22.9 cm x 22.9 cm) drip pan
- 2 front storage drawers
- 2 side storage drawers

Dimensions

- Overall: 27" W x 58.5" - 77.5" L x 31.75" H
(68.58cm W x 148.59 - 196.85cm L x 80.65cm H)
- Stirrup: 3 lateral positions 22°, 32°, 42° infinitely adjustable over length of bar
- Front drawers: (2) 18.25" W x 15.25" D x 4" H
(48.26cm W x 40.64cm D x 10.79cm H)
- Side drawers: (2) 18.25" W x 15.25" D x 7" H
(48.26cm W x 40.64cm D x 17.78cm H)



Fully-articulating back to position your patient as needed

Options

- 035 Articulated armboard assembly (specify patient R or L)
- 039 Accessory rail (seat or back section) (factory installed)
- 251 Bierhoff knee crutches (pair)
- 350 Pelvic tilt 7.5°
- 375 Safety side rails - pair (requires 2 each 039)
- 810 Restraint straps (set)
- 820 Stainless steel drain pan
- 851 Warmer drawer

Accent and upholstery colors

NOTE: Base color is only available in alabaster



United Metal Fabricators, Inc.

1316 Eisenhower Blvd. • Johnstown, PA 15904

800-638-5322 • 814-266-1870 fax • www.umf-exam.com